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
I. CONGRATULATIONS

Thank you for choosing Seinbiose to accompany you on your journey through motherhood! At Seinbiose we know that breast milk is the best natural food for your baby's healthy development and growth. Breastfeeding provides your baby with many natural benefits and allows you to enjoy intimate and tender moments together. As your trusted partner, Seinbiose helps mothers give the best to their babies where circumstances make breastfeeding impossible. Discover the ease and benefits of Seinbiose for you and your baby.

II. PRODUCT DESCRIPTION AND FEATURES

- Electric breast pump with easy-to-read LCD display.
- Double pumping mode.
- Two breast pumping phases: Stimulation and Expression.
- Stimulation phase : 5 levels / Expression phase : 9 levels.
- Compact and light-weight for maximum portability.
- Silent operation.
- Powered by AC adapter or 6xAAA batteries.

1- Dust protector	2 PP
2- Silicone breast shield	2 SIL
3- Plastic brace	2 PP
4- Diaphragm	2 SIL
5- Cap	2 PP
6- Cap plug	2 SIL
7- Bottle connector	2 PP
8- Tubing connector	4 PP/TPE
9- Tubing	2 SIL
10- Pump unit	ABS
11- AC Adapter	1 PP
12- Bottle lid	2 PP
13- Valve membrane	2 SIL
14- Valve	2 PP
15- Wide-neck milk bottle	2 PP
16- Bottle stand	2 PP

SIL = Silicone

PP = Polypropylene

ABS : Acrylonitrile butadiene styrene

III. IMPORTANT SAFETY INSTRUCTIONS

Please read through the safety instructions carefully before using your Seinbiose breast pump for the first time. Keep these safety instructions for future reference.

- Make sure that all parts are correctly assembled before use.
- Before use, clean and dry all parts in direct contact with the breast milk.
- Do not clean the cap unless necessary.
- Verify that there is no breast milk or other liquids in the tubing. Ensure that the tubing is completely dry before use.
- This device is for personal use. Use by multiple persons may present a health risk.
- Make sure that the voltage marked on the AC adapter is compatible with the power supply.
- Do not leave the breast pump unattended when plugged into a power source.

- Unplug the AC adapter from the power source when the pump is not in use, during cleaning and maintenance.
- Keep the breast pump unit and its parts out of the reach of babies and infants.
- Keep the breast pump unit and its parts away from heat sources and hot surfaces.
- Do not place the breast pump unit near appliances generating a strong magnetic field.
- Keep the breast pump unit dry and do not immerse in water. Do not spray or pour liquids over the breast pump.
- Do not touch an electrical device that has fallen into water. Unplug immediately.
- Do not use the breast pump in the bath or shower.
- Do not use the breast pump when sleepy or drowsy.
- Avoid exposing the breast pump and its parts to direct sunlight.
- Do not insert any foreign object into any of the openings of the breast pump and its accessories.
- Do not use the breast pump if the cable or pump unit is damaged.
- Do not dry or sterilise the breast pump or its parts in the microwave.
- Only use the AC power adapter provided by the manufacturer.
- Discontinue using your Seinbiose breast pump if you experience discomfort or pain during the pumping process and consult your healthcare professional.

IV. PREPARING THE PUMP


Ensure that all parts of the breast pump kit that are in direct contact with breast milk are cleaned or sterilised before each use: silicone breast shield (2), bottle connector (7), valve membrane (13), valve (14), milk bottle (15).

MOUNTING :

a. Assembling the breast pump kit

1. Positioning the silicone breast shield

Firmly press the silicone breast shield (2) into opening of the plastic brace (3) and attach to the bottle connector. Cover the silicone breast shield with the dust cover (1) when not in use.

2. Positioning the diaphragm

Place the diaphragm (4) into the opening at the top of the bottle connector (7). Cover the diaphragm with the cap (5). Seal tightly to avoid air leaks.

3. Assembling the valve membrane

Place the valve membrane (14) over the valve opening (13). Ensure that the membrane completely seals the valve opening to avoid air leaks. Attach the valve and valve membrane (14+13) to the bottom of the bottle connector (7) and seal tightly.

4. Assembling the bottle and bottle stand


Screw the bottle connector (7) onto the bottle (15) and place in the bottle stand (16).

5. Repeat steps 1-4 to assemble the second breast pump kit.

b. Connecting the breast pump kit to the breast pump unit

1. Attaching the breast pump kit to the pump unit

Connect one end of the tubing (9) to the breast pump kit and the other end to the pump unit (10), as illustrated. Repeat to connect the second breast pump kit. Attach the cable of the AC adapter (11) to the pump unit (10) and plug the AC adapter into a power outlet.

2. Using AA 1.5V batteries

Lift the cover of the battery compartment. Insert six (6) AA 1.5 V batteries in the direction indicated. Replace the cover.

V. USING YOUR SEINBIOSE BREAST PUMP


a. Before use:

- Wash your hands before handling the breast pump kit and unit.
- Verify that the breast pump kit is clean and dry.
- Verify that the tubing is dry. Do not operate your breast pump if the tubing contains water, liquids or other remnants. Clean and dry prior to use.
- Verify that your breast pump kit is assembled correctly and securely.

CAUTION: Only use the power adapter supplied with the electric breast pump. Make sure the power adapter is compatible with the power source.

If you are using your breast pump in battery mode, ensure that you are using appropriate batteries and that the batteries have been inserted correctly.

b. Understanding your display

c. Using your breast pump

- Find a calm and comfortable place that will allow you to draw your breast milk without being disturbed.
- Place the breast shield onto your breast. Your nipple should be at the centre of the breast shield tunnel to avoid rubbing against the sides. To ensure optimal suction, verify that your breast shield comfortably seals around your breast.
- Press the On/Off button to turn your breast pump on. The LCD screen will light up briefly. Press the on On/Off button again to start using the pump. The 30-minute integrated timer will begin immediately*.
- Your Seinbiose breast pump has two phase settings: Stimulation phase and Expression phase. The pump will automatically begin in Stimulation phase for 2 minutes before entering the Expression phase (factory default). You can also switch between the two phases manually by pressing the Phase button (). Apply the Stimulation phase until your breast milk begins to flow (approx. 2 minutes), then press the Phase button to switch to the Expression phase.
- To adjust the suction intensity to your comfort level during pumping, press the (+) to increase the suction, or the (-) to decrease the suction. Always begin with the most gentle and mild suction settings and then adjust to your needs and comfort level. Choosing the strongest setting immediately may result in lesions.
- Approximately 2 minutes after you cease to eject or express milk, stop and unplug the pump.

***Timing function:**

- The timer will automatically begin when you press the On/Off button of your breast pump. The timer is set to 30 minutes. Your breast pump will automatically turn off and enter standby mode once the 30-minute cycle has expired.

d. After use:

- Once finished, turn off your breast pump and unplug the device.
- Disassemble and clean your breast pump kit (all parts in direct contact with your breast milk). To clean, follow the cleaning instructions outlined in this user manual.

e. Storing and thawing breast milk

Maternal breast milk provides your baby with the right amount of vitamins, mineral salts, trace elements, sugars, fats and proteins in order to grow and develop in a healthy way. Store and thaw your breast milk appropriately in order to preserve your breast milk's nutritional qualities and to minimise bacterial contamination.

Breast milk can be preserved in the collection cup or bottle provided as part of your breast pump kit, or in freezer-safe containers and bags. Mark the date and hour at which you collected your breast milk on the container before storing it.

Seinbiose recommends the following guidelines for storing and thawing your breast milk :

Breast milk	Room temperature < 26 ° C	Refrigerator + 4 ° C	Freezer - 18 ° C
Freshly expressed milk	4 hours	48 hours	4 months
Thawed in the refrigerator but not heated	Max. 1 hour	24 hours	Do not re-freeze
Thawed and heated	Max. 1 hour	Max. 1 hour	Do not re-freeze
Milk remaining in the feeding bottle	To be consumed within 30 minutes	Throw out	Throw out

Tips for freezing and thawing breast milk:

- Never fill bottles / bags more than $\frac{3}{4}$ full to allow space for possible expansion.
- Label the bottle / bag with the date on which you expressed the milk.
- Do not store breast milk in the refrigerator door. Store it in the lower back area of the refrigerator, where the temperature is coldest.
- Do not heat your breast milk in the microwave. Microwaves may burn your breast milk or alter its nutritional qualities.
- Thaw frozen breast milk in the refrigerator overnight or place it under lukewarm running water (max. 37 degrees). Do not thaw frozen breast milk in boiling water.
- Shake the bottle or bag gently to blend any fat that has separated
- Use a small cooler with cool-packs to transport your breast milk. Do not exceed a maximal transportation duration of 1 hour.
- For premature or hospitalised babies, consult the healthcare professional accompanying your baby's development or your local milk bank for milk preservation and storage advice.

VI. HOW TO CLEAN YOUR BREAST PUMP

Before using your Seinbiose Double breast pump, disassemble and wash all parts of the breast pump kit that are in direct contact with the breast milk. Do not wash the tubing and the tubing connector. Completely air-dry all parts prior to next use.

IMPORTANT: Partially dry components may affect the operation of your electric breast pump. Ensure that all components are completely dry before use.

General cleaning instructions:

- Set aside the tubing and tubing connector. Under normal circumstances, the valve membrane will prevent milk backflow into the tubing, therefore these items do not need to be cleaned. However, should you find milk or water or other remnants in the tubing, follow the below Tubing Cleaning Instructions.
- All parts of your breast pump kit (**excluding the tubing and pump unit**) can be washed using a dishwasher or by hand using warm soapy water and a soft-bristle brush. Apply a complete washing cycle at 65 degrees Celsius.
- Do not use solvents or abrasives.
- Rinse with clean water.
- To clean the electric pump unit, wipe the unit with a damp cloth using only water or water mixed with a mild soap.

Tubing Cleaning Instructions:

- Wipe the outside of the tubing with a clean, damp cloth. Avoid getting water or other fluids inside the tubing.
- If you find liquids inside the tubing, rinse the tubing with clean hot water. Hang the tubing to air-dry completely before use. Do not insert any object into the tubing.

IMPORTANT: Prior to use, visually inspect the individual components of the breast pump kit for cracks, chips, tears, discoloration or general damage and deterioration. Damage to the various components may affect the efficiency of your Seinbiose Double breast pump. Discontinue use until all parts have been replaced.

VII. TROUBLESHOOTING

There is no or little suction

- Verify that all connections are properly attached at all connection points.
- Verify that the valve membrane is flat and completely seals the opening of the valve.
- Verify that the cap plug is correctly attached to the cap.

- Verify that there is no water in the tubing and that the tubing is not bent.
- Ensure that the silicone breast shield forms a complete seal around the breast.

The pump does not turn on

- Verify that the power chord is properly connected to the pump unit.
- Verify that the AC adapter is correctly inserted into a functional power source.
- Verify that the batteries work.
- Contact your dealer or the store where you bought your Seinbiose Double breast pump.

The pump is wet

- Unplug the AC adapter from the power source, then switch off the pump.
- Do not place the pump upside down; the control panel must always face upwards.
- Store the pump in a dry warm place for 24 hours

Note:

- Do not expose the pump to direct sunlight.
- Do not accelerate the drying process by use of heat.

Milk backflow or overflow

- If milk or other liquids flow into the tubing, stop your breast pump immediately. Unplug the pump unit and clean the tubing as described above. Leave the tubing to dry completely before use.
- Verify that the valve membrane is intact. If broken, replace it with a new one.

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